

Apc Mouse Anti-Human Cd45Ro Monoclonal Antibody For Fc - P08575

Item Number: CM0012040



Introduction

APC-conjugated mouse monoclonal antibody specific for human CD45RO, an isoform of the leukocyte common antigen PTPRC. Suited for flow cytometry (FC) in T-cell activation and immune profiling studies.

[Learn More](#)

Parameter	Value
Product Name	APC Mouse anti-Human CD45RO mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
Gene ID (Human)	5788
Gene ID	5788
Immunogen	Synthetic peptide A synthesized peptide derived from Human CD45RO
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08575
Protein Weight	136kDa
Shipping	Ice bag