

Apc Mouse Anti-Human Cd56 Monoclonal Antibody For Fc - P13591

Item Number: CM0012255



Introduction

APC-conjugated mouse monoclonal antibody targeting human CD56 (NCAM1), validated for flow cytometry. Suitable for studying neural cell adhesion, viral entry mechanisms, and NK cell biology. Reacts with human protein isoforms ranging from 40 to 95 kDa.

[Learn More](#)

Parameter	Value
Product Name	APC Mouse anti-Human CD56 mAb
Remarks/Alias	CD56; NCAM; MSK39
Species	Human
Gene ID	4684
Immunogen	Synthetic peptide derived from human CD56
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P13591
Protein Weight	40kDa/73kDa/80kDa/83kDa/93kDa/95kDa
Shipping	Ice bag