

Apc Mouse Anti-Human Cd62P/P-Selectin Mab - P16109

Item Number: CM0012323



Introduction

APC-conjugated mouse monoclonal antibody targeting human CD62P/P-selectin. Validated for flow cytometry, suitable for studying leukocyte rolling, cell adhesion, and inflammation in endothelial and platelet models.

[Learn More](#)

Parameter	Specification
Product Name	APC Mouse anti-Human CD62P/P-selectin mAb
Remarks/Alias	CD62; GRMP; PSEL; CD62P; GMP140; LECAM3; PADGEM
Species	Human
GeneID (Human)	6403
GeneID	6403
Immunogen	Recombinant protein of human CD62P/P-selectin.
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P16109
Protein Weight	91kDa
Shipping	Ice bag