

ApC Mouse Anti-Human Tcr A/B Monoclonal Antibody For Fc - P01848

Item Number: CM0004324



Introduction

APC-conjugated mouse monoclonal antibody targeting the constant region of human TCR alpha chain (TRAC). Validated for flow cytometry (FC) detection of human TCR α/β + T cells. Ideal for immunophenotyping and T cell receptor studies.

[Learn More](#)

Parameter	Specification
Product Name	APC Mouse anti-Human TCR α/β mAb
Remarks/Alias	IMD7; TCRA; TRA@
Species	Human
Gene ID (Human)	6955
Gene ID	6955
Immunogen	Recombinant protein of human TCR α/β
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P01848
Protein Weight	16kDa
Shipping	Ice bag