

Apc Rabbit Anti Ddddk-Tag Mab - Ddddk Tag

Item Number: CM0005907

Introduction

APC-conjugated rabbit monoclonal antibody against the DDDDK tag, suitable for flow cytometry (FC) including intracellular staining. Species-independent cross-reactivity ensures broad utility in tagged protein detection and research applications.

[Learn More](#)

Parameter	Details
Product Name	APC Rabbit anti DDDDK-Tag mAb
Remarks/Alias	DDDDK; DDDDK tag; DDDDK-tag
Species	Human
Immunogen	Synthetic peptide corresponding to DDDDK tag
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC, FC (intra)
Cross Reactivity	Species independent
Shipping	Ice bag