

Apc Rabbit Anti-Goat Igg (Fc) Fragment Monoclonal Antibody

Item Number: CM0012313

Introduction

APC-conjugated rabbit monoclonal antibody specific for goat IgG (Fc) fragment. Validated for flow cytometry (FC) applications. Enables detection of goat immunoglobulins in multiplexed immunoassays and diagnostic research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit Anti-Goat IgG (Fc) fragment mAb
Species	Goat
Immunogen	Recombinant protein of goat IgG (Fc) fragment
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Goat
Shipping	Ice bag