


Apc Rabbit Anti-Human Annexin A1/Anxa1 Mab For Fc - P04083

Item Number: CM0008806



Introduction

APC-conjugated rabbit monoclonal antibody targeting human Annexin A1/ANXA1. Designed for flow cytometry applications with human samples. Annexin A1 is a calcium-dependent phospholipid-binding protein involved in anti-inflammatory responses, T-cell differentiation, and phagocytosis.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human Annexin A1/ANXA1 mAb
Remarks/Alias	ANX1; LPC1
Species	Human
Gene ID (Human)	301
Gene ID	301
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2-346 of human Annexin A1/ANXA1 (NP_000691.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P04083
Protein Weight	39 kDa
Shipping	Ice bag