

Apc Rabbit Anti-Human Cd10 Monoclonal Antibody For Flow Cytometry - P08473

Item Number: CM0008036



Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD10 (CALLA/neprilysin) for flow cytometry. Recognizes the MME-encoded metalloprotease involved in opioid peptide and ANF degradation. Validated in human samples.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD10 mAb
Remarks/Alias	NEP; SFE; CD10; CALLA; CMT2T; SCA43
Species	Human
Gene ID (Human)	4311
Gene ID	4311
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 52-750 of human CD10 (NP_000893.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08473
Protein Weight	86kDa
Shipping	Ice bag