

ApC Rabbit Anti-Human Cd117/C-Kit Monoclonal Antibody For Fc - P10721

Item Number: CM0008159



Introduction

APC-conjugated rabbit monoclonal antibody against human CD117/c-Kit (P10721) for flow cytometry. This 110kDa receptor tyrosine kinase is essential for hematopoiesis, stem cell maintenance, and mast cell development. Suitable for IVD research applications.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD117/c-Kit mAb
Remarks/Alias	PBT; SCFR; C-Kit; CD117; MASTC
Species	Human
GeneID	3815
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-520 of human CD117/c-Kit (NP_000213.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P10721
Protein Weight	110kDa
Shipping	Ice bag