

Apc Rabbit Anti-Human Cd11A Monoclonal Antibody For Fc - P20701

Item Number: CM0004467



Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD11a (Integrin alpha-L, LFA-1) for flow cytometry. Suitable for studying leukocyte adhesion, migration, and immune synapse formation.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD11a mAb
Remarks/Alias	CD11A; LFA-1; LFA1A
Species	Human
Gene ID (Human)	3683
Gene ID	3683
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 26-340 of human CD11a (NP_002200.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P20701
Protein Weight	118kDa/129kDa/134kDa
Shipping	Ice bag