


APC Rabbit Anti-Human CD123 Monoclonal Antibody For Flow Cytometry - P26951

Item Number: CM0008182



Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD123 (IL-3R α) for flow cytometry (FC). Suitable for hematopoietic research; recognizes human CD123 with a ~43 kDa protein weight.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD123 mAb
Remarks/Alias	IL3R; CD123; IL3RX; IL3RY; IL3RAY; hIL-3Ra; IL-3R-alpha
Species	Human
Gene ID	3563
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-378 of human CD123 (NP_002174.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P26951
Protein Weight	43 kDa
Shipping	Ice bag