

Apc Rabbit Anti-Human Cd153 Monoclonal Antibody For Fc - P32971

Item Number: CM0012385

Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD153 (CD30L, TNFSF8). Recognizes human CD153, validated for flow cytometry. Suitable for T-cell proliferation and CD30/CD153 interaction research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD153 mAb
Remarks/Alias	CD153; CD30L; CD30LG; TNLG3A
Species	Human
GenelD (human)	944
GenelD	944
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 63-234 of human CD153 (NP_001235.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P32971
Protein Weight	26kDa
Shipping	Ice bag