

APC Rabbit Anti-Human Cd162/PSgl-1 Monoclonal Antibody For Flow Cytometry - Q14242

Item Number: CM0008579



Introduction

APC-conjugated rabbit monoclonal antibody against human CD162/PSGL-1 (P-selectin glycoprotein ligand 1), validated for flow cytometry. Suitable for studying leukocyte rolling, adhesion, and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD162/PSGL-1 mAb
Remarks/Alias	CLA; CD162; PSGL1; PSGL-1
Species	Human
Gene ID	6404
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 42-295 of Human CD162/PSGL-1 (NP_002997.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q14242
Shipping	Ice bag