

Apc Rabbit Anti-Human Cd185/Cxcr5 Monoclonal Antibody For Fc - P32302

Item Number: CM0010017

Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD185/CXCR5 (BLR1), a chemokine receptor involved in B-cell migration and lymphoma. Validated for flow cytometry. ~42 kDa membrane protein.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD185/CXCR5 mAb
Remarks/Alias	BLR1; CD185; MDR15
Species	Human
Gene ID (Human)	643
Gene ID	643
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-47 of Human CD185/CXCR5 (NP_033981.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P32302
Protein Weight	42kDa
Shipping	Ice bag