

ApC Rabbit Anti-Human Cd196/Ccr6 Monoclonal Antibody For Fc - P51684

Item Number: CM0005074



Introduction

APC-conjugated rabbit monoclonal antibody against human CD196/CCR6, validated for flow cytometry (FC). Ideal for studying chemokine receptor-mediated immune cell trafficking and inflammation.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD196/CCR6 mAb
Remarks/Alias	BN-1; DCR2; DRY6; CCR-6; CD196; CKRL3; GPR29; CKR-L3; CMKBR6; GPRCY4; STRL22; CC-CKR-6; C-C CKR-6
Species	Human
GeneID (human)	1235
GeneID	1235
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P51684
Protein Weight	42kDa
Shipping	Ice bag