


APC Rabbit Anti-Human Cd1B Monoclonal Antibody For Flow Cytometry - P29016

Item Number: CM0005125



Introduction

APC Rabbit anti-Human CD1B Monoclonal Antibody for FC. Specific for human, validated for flow cytometry, targeting CD1B (P29016), a T-cell surface glycoprotein that presents lipid antigens.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD1B mAb
Remarks/Alias	R1; CD1; CD1A
Species	Human
GeneID (Human)	910
GeneID	910
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P29016
Protein Weight	37kDa/30kDa
Shipping	Ice bag