


ApC Rabbit Anti-Human Cd1C Monoclonal Antibody For Flow Cytometry - P29017

Item Number: CM0010711



Introduction

APC-conjugated recombinant rabbit monoclonal antibody against human CD1c (BDCA1), validated for flow cytometry. Detects a lipid-antigen presenting glycoprotein on immune cells, suitable for immunology and dendritic cell research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD1c mAb
Remarks/Alias	R7; CD1; CD1A; BDCA1
Species	Human
Gene ID (Human)	911
Gene ID	911
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-298 of human CD1c (NP_001756.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P29017
Protein Weight	38kDa
Shipping	Ice bag