


ApC Rabbit Anti-Human Cd209/Dc-Sign Monoclonal Antibody For Fc - Q9Nnx6

Item Number: CM0009803



Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD209/DC-SIGN (Q9NNX6). Validated for flow cytometry, it enables precise detection of dendritic cells and study of immune response, pathogen recognition, and host-pathogen interactions.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD209/DC-SIGN mAb
Remarks/Alias	CDSIGN; CLEC4L; DC-SIGN; DC-SIGN1; hDC-SIGN
Species	Human
Gene ID (Human)	30835
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 59-404 of human CD209/DC-SIGN (NP_066978.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9NNX6
Protein Weight	4kDa/18kDa/30kDa/33kDa/35kDa/37kDa/41kDa/43kDa/45kDa
Shipping	Ice bag