


Apc Rabbit Anti-Human Cd213A1/IL-13R α 1 Monoclonal Antibody For Fc - P78552

Item Number: CM0011514



Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD213a1/IL-13R α 1, validated for flow cytometry. Ideal for studying IL-13/IL-4 signaling, cancer/testis antigen expression, and cytokine receptor biology.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD213a1/IL-13R α 1 Monoclonal Antibody
Remarks/Alias	NR4; CT19; CD213A1; IL-13Ra
Species	Human
GeneID	3597
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-343 of human CD213a1/IL-13R α 1 (NP_001551.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P78552
Protein Weight	49kDa
Shipping	Ice bag