

Apc Rabbit Anti-Human Cd228/Mfi2 Monoclonal Antibody For Flow Cytometry - P08582

Item Number: CM0005482



Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD228/Melanotransferrin (MFI2). Validated for flow cytometry, suitable for research on iron uptake, melanoma, and related areas.

[Learn More](#)

Parameter	Specification
Product Name	APC Rabbit anti-Human CD228/MFI2 mAb
Remarks/Alias	MTf; MFI2; MTF1; CD228; MAP97
Species	Human
GeneID (human)	4241
GeneID	4241
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-708 of human CD228/MFI2 (NP_005920.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08582
Protein Weight	32kDa/80kDa
Shipping	Ice bag