

Apc Rabbit Anti-Human Cd244/2B4 Monoclonal Antibody For Flow Cytometry - Q9Bzw8

Item Number: CM0012800

Introduction

APC-conjugated rabbit monoclonal antibody against human CD244/2B4 (SLAMF4) for flow cytometry. Recognizes natural killer cell receptor 2B4 involved in NK cell activation. Validated in human samples.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD244/2B4 mAb
Remarks/Alias	2B4; NAIL; Nmrk; NKR2B4; SLAMF4
Species	Human
Gene ID (Human)	51744
Gene ID	51744
Immunogen	Recombinant protein corresponding to a sequence within amino acids 22-221 of human CD244/2B4 (NP_001160135.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9BZW8
Protein Weight	42kDa
Shipping	Ice bag