

APC Rabbit Anti-Human Cd261/Trail-R1 Monoclonal Antibody For Fc - 000220

Item Number: CM0012066



Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD261/TRAIL-R1 (DR4, TNFRSF10A) for flow cytometry. Recognizes a 50kDa death receptor involved in TRAIL-induced apoptosis and NF-kB activation. Validated for human reactivity.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD261/TRAIL-R1 mAb
Remarks/Alias	TNFRSF10A; APO2; CD261; DR4; TRAILR-1; TRAILR1; TNF receptor superfamily member 10a
Species	Human
GeneID (Human)	8797
GeneID	8797
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 124-232 of human CD261/TRAIL-R1 (NP_003835.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O00220
Protein Weight	50kDa
Shipping	Ice bag