

Apc Rabbit Anti-Human Cd27 Monoclonal Antibody For Fc - P26842

Item Number: CM0010999



Introduction

High-quality APC-conjugated rabbit anti-human CD27 monoclonal antibody (clone unspecified) validated for flow cytometry (FC). Suitable for studying costimulatory immune-checkpoint receptor expression on T-cells, NK-cells, and B-cells in immunology and oncology research.

[Learn More](#)

Parameter	Specification
Product Name	APC Rabbit anti-Human CD27 mAb
Remarks/Alias	T14; S152; Tp55; TNFRSF7; S152, LPFS2
Species	Human
GeneID	939
Immunogen	Recombinant fusion protein containing amino acids 21-192 of human CD27 (NP_001233.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P26842
Protein Weight	29kDa
Shipping	Ice bag