

## ApC Rabbit Anti-Human Cd270/Hvem Monoclonal Antibody For Flow Cytometry - Q92956

Item Number: CM0011155



### Introduction

APC-conjugated rabbit monoclonal antibody against human CD270/HVEM (TNFRSF14), designed for flow cytometry applications in immunology and virology research. Detects endogenous CD270 on immune cells.

[Learn More](#)

| Specification    | Details                                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | APC Rabbit anti-Human CD270/HVEM mAb                                                                                                        |
| Alias            | TNFRSF14; ATAR; CD270; HVEA; HVEM; LIGHTR; TR2; TNF receptor superfamily member 14                                                          |
| Species          | Human                                                                                                                                       |
| Gene ID (Human)  | 8764                                                                                                                                        |
| Immunogen        | Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 39-162 of human CD270/HVEM (NP_003811.2) |
| Source           | Rabbit                                                                                                                                      |
| Category         | Monoclonal Antibody                                                                                                                         |
| Application      | Flow Cytometry (FC)                                                                                                                         |
| Cross Reactivity | Human                                                                                                                                       |
| SWISS            | Q92956                                                                                                                                      |
| Protein Weight   | 21kDa/30kDa                                                                                                                                 |
| Shipping         | Ice bag                                                                                                                                     |