

Apc Rabbit Anti-Human Cd271/Ngfr Monoclonal Antibody For Fc - P08138

Item Number: CM0008271



Introduction

High-quality APC-conjugated rabbit monoclonal antibody against human CD271/NGFR (p75NTR) for flow cytometry applications. Ideal for studying neurotrophin signaling, neuronal development, and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD271/NGFR mAb
Remarks/Alias	CD271; p75NTR; TNFRSF16; p75(NTR); Gp80-LNGFR
Species	Human
GeneID (Human)	4804
GeneID	4804
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08138
Protein Weight	35kDa/45kDa
Shipping	Ice bag