

APC Rabbit Anti-Human Cd33 Monoclonal Antibody For Flow Cytometry - P20138

Item Number: CM0009523



Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD33 (Siglec-3, gene ID 945), validated for flow cytometry. Ideal for studying myeloid lineage cells, microglia, and AML blasts. SWISS: P20138.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD33 mAb
Remarks/Alias	p67; SIGLEC3; SIGLEC-3
Species	Human
GeneID	945
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-259 of human CD33 (NP_001763.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (FC)
Cross Reactivity	Human
SWISS	P20138
Protein Weight	40 kDa
Shipping	Ice bag