

Apc Rabbit Anti-Human Cd38 Monoclonal Antibody For Fc - P28907

Item Number: CM0012050



Introduction

APC-conjugated rabbit monoclonal antibody against human CD38, suitable for flow cytometry. Targets the ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 1 (P28907), involved in calcium signaling and lymphocyte activation.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD38 mAb
Remarks/Alias	ADPRC1; cADPR1; ADPRC 1
Species	Human
Gene ID (Human)	952
Gene ID	952
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 43-300 of human CD38 (P28907).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P28907
Protein Weight	13kDa/34kDa
Shipping	Ice bag