


ApC Rabbit Anti-Human Cd40 Monoclonal Antibody For Fc - P25942

Item Number: CM0012741



Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD40 (TNFRSF5, Bp50) validated for flow cytometry. Ideal for immunology, oncology research, and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD40 mAb
Alias	p50; Bp50; CDW40; TNFRSF5
Species	Human
Gene ID	958
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-193 of human CD40 (NP_001241.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross-Reactivity	Human
UniProt Accession	P25942
Calculated Molecular Weight	31kDa
Shipping Condition	Ice bag