

Apc Rabbit Anti-Human Cd42B Monoclonal Antibody For Fc - P07359

Item Number: CM0004068



Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD42b (GP1BA). Validated for flow cytometry (FC). Suitable for platelet function, thrombosis, and hemostasis research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD42b mAb
Remarks/Alias	GP1BA; BDPLT1; BDPLT3; BSS; CD42B; CD42b-alpha; DBPLT3; GP1B; GPIbA; GPIbalpha; VWDP; platelet glycoprotein Ib alpha chain
Species	Human
GeneID (Human)	2811
GeneID	2811
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 15-531 of human CD42b (NP_000164.5).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P07359
Protein Weight	72kDa
Shipping	Ice bag