


Apc Rabbit Anti-Human Cd44 Monoclonal Antibody For Fc - P16070

Item Number: CM0011986



Introduction

APC-conjugated recombinant rabbit monoclonal antibody targeting human CD44 (P16070), validated for flow cytometry. Ideal for lymphocyte homing, adhesion, and migration studies.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD44 mAb
Remarks/Alias	IN; LHR; MC56; MDU2; MDU3; MIC4; Pgp1; CDW44; CSPG8; H-CAM; HCELL; ECM-III; HUTCH-1; HUTCH-I; ECMR-III; Hermes-1
Species	Human
GeneID (Human)	960
GeneID	960
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-178 of human CD44 (NP_000601.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P16070
Protein Weight	82kDa
Shipping	Ice bag