

Apc Rabbit Anti-Human Cd45Rb Monoclonal Antibody For Fc - P08575-9

Item Number: CM0011049

Introduction

APC-conjugated rabbit monoclonal antibody against human CD45RB (PTPRC isoform 7) for flow cytometry. Detects human CD45RB, the leukocyte common antigen. Ideal for immune cell immunophenotyping research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD45RB mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
Gene ID (Human)	5788
Gene ID	5788
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 26-463 of human CD45RB (XP_006711537.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08575-9
Shipping	Ice bag