


Apc Rabbit Anti-Human Cd47 Monoclonal Antibody For Fc - Q08722

Item Number: CM0006242



Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD47 (IAP/OA3/MER6) for flow cytometry. Validated with recombinant immunogen (aa 19-139), reacts with human. Molecular weight 31-35kDa.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD47 mAb
Remarks/Alias	IAP; OA3; MER6
Species	Human
GeneID (human)	961
GeneID	961
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 19-139 of human CD47 (NP_001768.1))
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q08722
Protein Weight	31kDa/33kDa/35kDa
Shipping	Ice bag