

APC Rabbit Anti-Human Cd48 Monoclonal Antibody, Glycerol-Free For Flow Cytometry - P09326

Item Number: CM0009349



Introduction

APC-conjugated rabbit monoclonal antibody against human CD48 (SLAMF2), glycerol-free, validated for flow cytometry. Suitable for immunophenotyping of hematopoietic cells and functional studies of CD48-mediated immune interactions.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD48 mAb, glycerol free
Remarks/Alias	BCM1; BLAST; hCD48; mCD48; BLAST1; SLAMF2; MEM-102
Species	Human
GeneID (Human)	962
GeneID	962
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-220 of human CD48(NP_001769.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P09326
Protein Weight	28kDa
Shipping	Ice bag