

ApC Rabbit Anti-Human Cd66A/Ceacam1 Monoclonal Antibody For Flow Cytometry - P13688

Item Number: CM0010968



Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD66a/CEACAM1, validated for flow cytometry. Useful for immune cell profiling, T cell, NK cell, and tumor cell studies. Recognizes human CEACAM1.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD66a/CEACAM1 mAb
Remarks/Alias	CEACAM1; BGP; BGP1; BGPI; carcinoembryonic antigen-related cell adhesion molecule 1
Species	Human
Gene ID (Human)	634
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-428 of human CD66a/CEACAM1 (NP_001703.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P13688
Protein Weight	58kDa
Shipping	Ice bag