


ApC Rabbit Anti-Human Cd81 Monoclonal Antibody For Flow Cytometry - P60033

Item Number: CM0011598



Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD81 (TAPA-1). Validated for flow cytometry, this antibody detects a tetraspanin involved in B-cell signaling and viral entry. Ideal for immunology and virology research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD81 mAb
Remarks/Alias	S5.7; CVID6; TAPA1; TSPAN28
Species	Human
GeneID (Human)	975
GeneID	975
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 114-199 of human CD81 (NP_004347.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P60033
Protein Weight	26kDa
Shipping	Ice bag