


Apc Rabbit Anti-Human Cd9 Monoclonal Antibody For Fc - P21926

Item Number: CM0008442



Introduction

APC-conjugated rabbit monoclonal antibody against human CD9, a tetraspanin involved in cell adhesion, platelet activation, and tumor metastasis. Validated for flow cytometry (FC) applications.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD9 mAb
Remarks/Alias	MIC3; MRP-1; BTCC-1; DRAP-27; TSPAN29; TSPAN-29
Species	Human
GeneID (Human)	928
GeneID	928
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-228 of human CD9 (NP_001760.1))
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P21926
Protein Weight	25kDa
Shipping	Ice bag