

Apc Rabbit Anti-Human Cd95/Fas Monoclonal Antibody For Flow Cytometry - P25445

Item Number: CM0005956



Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD95/FAS (Apo-1, TNFRSF6). Validated for flow cytometry, suitable for apoptosis studies, peripheral tolerance, and T-cell homeostasis research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD95/FAS mAb
Remarks/Alias	APT1; CD95; FAS1; APO-1; FASTM; ALPS1A; TNFRSF6
Species	Human
GeneID (Human)	355
GeneID	355
Immunogen	Recombinant Protein (Recombinant fusion protein containing a sequence corresponding to amino acids 26-173 of human CD95/FAS (NP_000034.1))
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P25445
Protein Weight	38kDa
Shipping	Ice bag