

ApC Rabbit Anti-Human Egfr Monoclonal Antibody For Flow Cytometry - P00533

Item Number: CM0011294



Introduction

APC-conjugated rabbit monoclonal antibody targeting human EGFR, validated for flow cytometry. Recombinant immunogen (aa 335–525). Ideal for oncology and signaling research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human EGFR mAb
Remarks/Alias	ERBB; ERBP; HER1; mENA; ERBB1; PIG61; NISBD2
Species	Human
GeneID (Human)	1956
GeneID	1956
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 335-525 of human EGFR (NP_005219.2))
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P00533
Protein Weight	44kDa/69kDa/77kDa/134kDa
Shipping	Ice bag