


Apc Rabbit Anti-Human Flt3/Cd135 Monoclonal Antibody For Fc - P36888

Item Number: CM0005442



Introduction

High-quality APC-conjugated rabbit monoclonal antibody against human Flt3/CD135 (FLK2/STK1). Validated for flow cytometry. Targets extracellular domain epitope; ideal for hematopoietic and leukemia research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human Flt3/CD135 Monoclonal Antibody
Remarks/Alias	FLK2; STK1; CD135; FLK-2
Species	Human
GeneID (Human)	2322
GeneID	2322
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-163 of human Flt3/CD135 (NP_004110.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P36888
Protein Weight	108KDa/113kDa
Shipping	Ice bag