


Apc Rabbit Anti-Human Icam-1/Cd54 Monoclonal Antibody For Flow Cytometry (Fc) - P05362

Item Number: CM0008624



Introduction

APC-conjugated rabbit monoclonal antibody targeting human ICAM-1 (CD54), validated for flow cytometry. Recombinant immunogen spanning amino acids 28-480. Suitable for studying leukocyte adhesion, transendothelial migration, and viral receptor function.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human ICAM-1/CD54 mAb
Remarks/Alias	BB2; CD54; P3.58
Species	Human
GeneID (Human)	3383
GeneID	3383
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-480 of human ICAM-1/CD54 (NP_000192.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P05362
Protein Weight	58kDa
Shipping	Ice bag