


ApC Rabbit Anti-Human IgG (Fc) Monoclonal Antibody For Flow Cytometry - P01857

Item Number: CM0009241



Introduction

APC-conjugated rabbit monoclonal antibody targeting human IgG Fc region, validated for flow cytometry. Ideal for B-cell research, humoral immunity studies, and therapeutic antibody development. - P01857

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human IgG (Fc) mAb
Remarks/Alias	IGHG1
Species	Human
Gene ID (Human)	3500
Gene ID	3500
Immunogen	Recombinant protein (or fragment). This information is considered to be commercially sensitive.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P01857
Protein Weight	36 kDa
Shipping	Ice bag