


ApC Rabbit Anti-Human IL-17A Monoclonal Antibody For Flow Cytometry (Intracellular) - Q16552

Item Number: CM0008202



Introduction

APC-conjugated rabbit monoclonal antibody against human IL-17A, validated for intracellular flow cytometry. Suitable for human samples; ideal for Th17 cell, inflammation, and mucosal immunity research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human IL-17A mAb
Remarks/Alias	IL17; CTLA8; IL-17; ILA17; CTLA-8; IL-17A
Species	Human
GeneID (Human)	3605
GeneID	3605
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-155 of human IL-17A (NP_002181.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	Q16552
Protein Weight	18kDa
Shipping	Ice bag