

Apc Rabbit Anti-Human IL-2 Mab - P60568

Item Number: CM0008187



Introduction

APC-conjugated rabbit monoclonal antibody against human IL-2 (Interleukin-2, TCGF). Designed for intracellular flow cytometry (FC). Recognizes human IL-2, essential for T-cell growth, immune regulation, and NK cell activity.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human IL-2 mAb
Remarks/Alias	IL-2; TCGF; lymphokine
Species	Human
GeneID (Human)	3558
GeneID	3558
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 21-153 of human IL-2 (NP_000577.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P60568
Protein Weight	18kDa
Shipping	Ice bag