

APC Rabbit Anti-Human Integrin-β1/CD29 Monoclonal Antibody For Flow Cytometry - P05556

Item Number: CM0006258



Introduction

APC-conjugated rabbit monoclonal antibody targeting human Integrin-β1/CD29 (UniProt P05556), validated for flow cytometry. Suitable for studying cell adhesion, migration, and receptor-mediated signaling in research applications.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human Integrin-β1/CD29 mAb
Remarks/Alias	CD29; FNRB; MDF2; VLAB; GPIIA; MSK12; VLA-BETA
Species	Human
GeneID (Human)	3688
GeneID	3688
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-728 of human Integrin-β1/CD29 (NP_002202.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P05556
Protein Weight	87kDa/88kDa/91kDa
Shipping	Ice bag