


Apc Rabbit Anti-Human Mertk Monoclonal Antibody For Fc - Q12866

Item Number: CM0004944



Introduction

APC-conjugated rabbit monoclonal antibody against human MERTK receptor tyrosine kinase. Validated for flow cytometry (FC) applications; cross-reacts with human. Suitable for studying efferocytosis, macrophage clearance, and cancer signaling.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human MERTK mAb
Remarks/Alias	MERTK; MER; RP38; Tyro12; c-Eyk; c-mer; tyrosine-protein kinase Mer
Species	Human
GeneID (Human)	10461
GeneID	10461
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-499 of human MERTK (NP_006334.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q12866
Protein Weight	110kDa
Shipping	Ice bag