

Apc Rabbit Anti-Human Mglur2 Monoclonal Antibody For Fc - Q14416

Item Number: CM0004586



Introduction

APC-conjugated rabbit monoclonal antibody against human metabotropic glutamate receptor 2 (mGluR2, GRM2) for flow cytometry. Validated for human samples. Suitable for neuroscience research on synaptic transmission and viral entry studies.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human mGluR2 mAb
Remarks/Alias	GLUR2; mGlu2; GPRC1B; MGLUR2
Species	Human
Gene ID (Human)	2912
Gene ID	2912
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 22-498 of human mGluR2 (NP_000830.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q14416
Protein Weight	96kDa
Shipping	Ice bag