

APC Rabbit Anti-Human Nectin-1/Cd111 Monoclonal Antibody For Flow Cytometry - Q15223

Item Number: CM0009599



Introduction

APC-conjugated rabbit monoclonal antibody targeting human Nectin-1/CD111. Validated for flow cytometry, ideal for cell adhesion, neurodevelopment, and viral entry studies. Reacts with human samples. Supports research in neuroscience, virology, and cell biology.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human Nectin-1/CD111 mAb
Remarks/Alias	ED4; PRR; HlgR; HV1S; HVEC; OFC7; PRR1; PVRR; CD111; PVRL1; PVRR1; SK-12; CLPED1; nectin-1
Species	Human
GeneID (Human)	5818
GeneID	5818
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-334 of human Nectin-1/CD111 (NP_002846.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q15223
Protein Weight	57kDa/51kDa/39kDa
Shipping	Ice bag