


Apc Rabbit Anti-Human Podoplanin Monoclonal Antibody For Flow Cytometry - Q86YI7

Item Number: CM0005802



Introduction

High-quality APC-conjugated rabbit monoclonal antibody targeting human podoplanin (PDPN), validated for flow cytometry applications. Essential for investigating cell adhesion, migration, and tumor invasion mechanisms. Supplied by CamelBio for diagnostic and research use.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human Podoplanin mAb
Remarks/Alias	T1A; GP36; GP40; Gp38; OTS8; T1A2; T11A; D2-40; T1A-2; AGGRUS; HT1A-1; PA2.26
Species	Human
GeneID (human)	10630
GeneID	10630
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-123 of Human Podoplanin (NP_006465.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q86YL7
Protein Weight	17kDa
Shipping	Ice bag