

ApC Rabbit Anti-Human Siglec-8 Monoclonal Antibody For Fc - Q9Nyz4

Item Number: CM0009310



Introduction

APC-conjugated rabbit monoclonal antibody targeting human Siglec-8, validated for flow cytometry. Recognizes sialic acid-binding lectin 8 on basophils, mast cells, and eosinophils. Research use only.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human Siglec-8 mAb
Remarks/Alias	SAF2; SIGLEC-8; SIGLEC8L
Species	Human
Gene ID (Human)	27181
Gene ID	27181
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-363 of human Siglec-8 (NP_055257.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9NYZ4
Protein Weight	44kDa/46kDa/54kDa
Shipping	Ice bag