


ApC Rabbit Anti-Human/Monkey Cd14 Monoclonal Antibody For Flow Cytometry - P08571

Item Number: CM0009938



Introduction

APC-conjugated rabbit monoclonal antibody recognizes human and cynomolgus CD14, suitable for flow cytometry. CD14 is a key innate immunity coreceptor for bacterial LPS, expressed on monocytes and macrophages.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human/Monkey CD14 mAb
Remarks/Alias	CD14
Species	Human
Gene ID (Human)	929
Gene ID	929
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-352 of human CD14 (NP_000582.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08571
Protein Weight	40kDa
Shipping	Ice bag